


NORVI

# TELI Ai

Automated Fever Screening with AI



Temperature Reading



Automated Alerts

33.7 °C

37.6 °C






34.2 °C

36.1 °C

34.4 °C

33.9 °C


33.1 °C



Product Brochure



Automated Alerts



Temperature Reading



## Features

Multi-person detection simultaneously

Detects faces with masks

Visible image of the person

Immediately alarming to notify operators

Zero physical contact measurement

## Applications

Airport

Coach Station

Railway

Freeway Exit

Mall Entrance

Bus



21/03/2020 14:32:11  
High Fever  
35.6 °c



# Models

## TELI Ai

AT12-1219TH

### High Accuracy Fever monitoring

Thermal resolution 160 x 120

Ethernet Connectivity

Measurement Accuracy 0.3 'c

2 x Alarm Outputs

Real-time alerts

micro-SD card storage

AT12-8060TH

### Cost effective fever monitoring

Thermal resolution 80 x 60

Ethernet / WiFi Connectivity

Measurement Accuracy 0.5 'c

2 x Alarm Outputs

Real-time alerts

E-mail alerts / Reporting

For further inquiries and tech support

**info@icd.lk**

ICONIC DEVICES PRIVATE LIMITED

Head office  
183, Maharagama Road  
Boralasgamuwa, Sri Lanka

The information provided in this documentation contains general descriptions and/or technical characteristics of the performance of the products contained herein. This documentation is not intended as a substitute for and is not to be used for determining suitability or reliability of these products for specific user applications. It is the duty of any such user or integrator to perform the appropriate and complete risk analysis, evaluation and testing of the products with respect to the relevant specific application or use thereof. Neither ICONIC DEVICES nor any of its affiliates or subsidiaries shall be responsible or liable for misuse of the information contained herein.